


 **DIAGMED**

HEALTHCARE



## Mouthguards

The widest choice of both disposable and re-usable mouthguards available in the UK





Diagmed Healthcare Ltd is a specialist distributor of world leading medical technology.

Diagmed's mission is to capitalise on our experience and assets by identifying and representing unique and innovative product lines that benefit patients and medical staff.

We fully understand exactly how to distribute medical/surgical devices within the UK health care markets. We share a good working relationship with a very large customer base and fully support our customers before and after the point of sale.

Our team of sales representatives cover the whole of the U.K, and undergo regular training on product knowledge and have the full technical support of our team of administrators at our head office.

## Standard Disposable Mouthguard

Features a large port with large lateral ports to aid suctioning and finger assisted scope insertion.

Product Code	Description
MG-1	Disposable Mouthguard
MG-1/S	Disposable Mouthguard with latex-free strap



## Oxyguard Disposable and Reusable Mouthguards



**Disposable.** Unique oxygen nasal port, large insertion port, comfortable, patient friendly.

**Reusable.** Unique oxygen nasal port, large insertion port, comfortable, autoclavable.

Product Code	Description
Oxyguard-D	Disposable Mouthguard
Oxyguard-DS	Disposable Mouthguard with Strap
Oxyguard	Reusable Mouthguard



## Disposable Mouthguards with Dental Retention Rim

Comfortable large port disposable mouthguards and coordinating straps are latex-free, addressing the ever-growing problem of latex sensitivity and risk of severe allergic reaction. This extensive line includes a variety of sizes, designs and models to ensure the best fit for the patient and clinician. An impressive selection – latex-free and worry free.

### Standard Clearview Mouthguard

Comfortable large port disposable mouthguard.

Product Code	Description
00712801	Standard Mouthguard with latex-free rubber strap
00712802	Standard Mouthguard without strap
00712803	Standard Mouthguard with velcro strap



### Maxi Mouthguard

Product Code	Description
00712804	Maxi Mouthguard with latex-free rubber strap
00712805	Maxi Mouthguard without strap
00712806	Maxi Mouthguard with velcro strap





## Disposable Mouthguards without Dental Retention Rim

Comfortable large port disposable mouthguards and coordinating straps are latex-free, addressing the ever-growing problem of latex sensitivity and risk of severe allergic reaction. This extensive line includes a variety of sizes, designs and models to ensure the best fit for the patient and clinician. An impressive selection – latex-free and worry free.



Product Code	Description
00712808	Paediatric Mouthguard with latex-free rubber strap
00712809	Intermediate Mouthguard with latex-free rubber strap



## G.I.S Bloc Mouthguards



### OxyBloc

OxyBloc combines the UltimaBloc with a nasal cannula and seven feet of oxygen tubing.

OxyBloc comes with a clip that attaches to the patient's garment and a latex-free headstrap to help keep the biteblock securely in place.



### UltimaBloc

Featuring a specially engineered Slip-Slope design, UltimaBloc naturally slides back into the mouth as the patient bites down.

The large lateral ports enable easy suctioning and finger assisted scope insertion. UltimaBloc accepts up to a 50 French dilator. Latex-free headstrap is included.

Product Code	Description
GIS-10D	OxyBloc Mouthguard
GIS-4D	UltimaBloc Mouthguard



## GeneraBloc

GeneraBloc improves on the traditional design of a disposable biteblock with a larger port.

This low-cost biteblock features a smooth port that accepts up to a 54 French dilator.



## OmniBloc

OmniBloc features a large port and a latex-free headstrap at a low cost.

Its high front flange maximizes patient comfort by keeping the biteblock from getting under the patient's lips. OmniBloc easily accepts a 60 French dilator.



## PediaBloc

PediaBloc has been designed specifically for use with small children, infants and the elderly.

It features a large flange that prevents the biteblock from slipping into the patient's mouth. Latex-free headstrap is included.



Product Code	Description
GIS-3D	GeneraBloc Mouthguard
GIS-11M	OmniBloc Mouthguard
GIS-2D	PediaBloc Mouthguard



PROGRESS AND EFFICIENCY THROUGH INNOVATION



**Diagmed Healthcare Ltd**

Lumley Close, Thirsk Industrial Park, Thirsk, North Yorkshire, England YO7 3TD

**Tel:** (01845) 526660 **Fax:** (01845) 522199 **email:** [info@diagmed.co.uk](mailto:info@diagmed.co.uk) **www:** [www.diagmed.co.uk](http://www.diagmed.co.uk)