




Disposable Polypectomy Snares

We believe our wide choice of Polypectomy accessories are unmatched by any other company.



Diagmed Healthcare Ltd is a specialist distributor of world leading medical technology.

Diagmed's mission is to capitalise on our experience and assets by identifying and representing unique and innovative product lines that benefit patients and medical staff.

We fully understand exactly how to distribute medical/surgical devices within the UK health care markets. We share a good working relationship with a very large customer base and fully support our customers before and after the point of sale.



Our team of sales representatives cover the whole of the U.K, and undergo regular training on product knowledge and have the full technical support of our team of administrators at our head office.

Disposable polypectomy snares




PolypMaster

Diagmed Healthcare introduces the first complete line of disposable snares. Quality performance disposable snares in various loop sizes and shape configuration - a snare to meet any requirement.

Product Code	Description
D4401	Sub Micro Oval - 1.0cm Loop - 2.3mm x 240cm (short throw handle)
D4402	Standard Oval - 3.0cm Loop - 2.3mm x 240cm
D4404	Mini Oval - 2.5cm Loop - 2.3mm x 240cm (short throw handle)
D4406	Micro Oval - 2.0cm Loop - 2.3mm x 240cm (short throw handle)
D4408	Standard Hexagonal - 3.0cm Loop - 2.3mm x 240cm
D4410	Standard Crescent - 2.5cm Loop - 2.3mm x 240cm
D4411	Jumbo Oval - 3.5cm Loop - 2.3mm x 240cm
D4412	Mini Crescent - 1.5cm Loop - 2.3mm x 240cm
D4413	Jumbo Hexagonal - 3.5cm Loop - 2.3mm x 240cm
D4414	Micro Crescent - 1.0cm Loop - 2.3mm x 240cm
D4416	Modified Crescent - 3.0cm Loop - 2.3mm x 240cm
D4422	Paediatric Snare - 2.5cm Loop - 1.8mm x 180cm (short throw handle)



Olympus diathermy compatible.



Disposable polypectomy snares

Kudo Snare

- Symmetrical polypectomy snare.
- Monofilament wire.
- Olympus Diathermy Compatible.
- Controllable, variable opening width

Product Code	Description
99-05-20-11-23-8	Standard Oval 40mm



Shark Tooth Snare

- Monofilament wire.
- 4 spikes for improved hold on tissue.
- Olympus Diathermy Compatible.

Product Code	Description
99-05-20-51-22-5	Standard Oval 25mm



Endo-Snare

- Short throw handle with centimetre markings for measuring polyp size.
- Olympus Diathermy Compatible.


Product Code	Description
SGH-1622C	Standard Hexagonal Snare 30mm
SGH-1622L	Standard Hexagonal Snare 30mm
SGO-1622S	Disposable Colonic Oval Hot Snare 10mm
SGO-1622M	Disposable Colonic Oval Hot Snare 10mm
SGO-1622L	Disposable Colonic Oval Hot Snare 10mm
SGO-1622XL	Disposable Colonic Oval Hot Snare 10mm



Disposable polypectomy snares

Carr-Locke and Grabber Snare

- A unique device for a specific procedure.
- Four unique stainless steel grabbers to secure (Sessile) polyp without tearing the mucosa.
- Single use eliminates washing/sterilizing.
- No time lost finding and assembling reusable equivalent snares.



Product Code	Description	Loop	Diameter Catheter	Length
D4432	Standard Carr-Locke Snare	3.0cm	2.3mm	240cm
D4450	Disposable Grabber Snare	2.5cm	2.3mm	240cm
D4454	Carr-Locke Grabber Snare	3.0cm	2.3mm	240cm

ProSnare

Polypectomy Snares with Braided and Monofilament Wire


Oval, Braided Wire

Product Code	Loop Size [Diameter x Length]	Diameter	Length
99-05-20-12-22-2 TW	10 x 20mm	2.3mm	230cm
99-05-20-12-20-2 TW	15 x 30mm	2.3mm	230cm
99-05-20-12-21-2 TW	20 x 40mm	2.3mm	230cm
99-05-20-12-24-2 TW	25 x 50mm	2.3mm	230cm
99-05-20-12-23-2 TW	30 x 55mm	2.3mm	230cm



Hexagonal, Braided Wire

Product Code	Loop Size [Diameter x Length]	Diameter	Length
99-05-20-12-23-4 TW	30 x 50mm	2.3mm	230cm



The following Diamed snares are made using monofilament wire, a stiffer wire improving ease of use when grasping mucosa.

Oval, Monofilament Wire

Product Code	Loop Size [Diameter x Length]	Diameter	Length
99-05-20-12-22-5 MW	20 x 45mm	2.3mm	230cm
99-05-20-12-23-5 MW	25 x 50mm	2.3mm	230cm




Braided WIRE

Disposable polypectomy snares

Skorpio Snare

The Diagmed Skorpio Snare combines injection needle and polypectomy snare for faster more efficient procedures.

- Inject or snare without device exchange
- Maintain scope position during polypectomy
- Monofilament snare position provides increased wire strength




Product Code	Opening Width (Diameter x Length)	Diameter	Length	Needle Diameter
99-05-21-11-22-5	20 x 40mm	2.6mm	230cm	0.7mm

Vario Snare

The unique Diagmed Vario Snare provides the user with standard and small diameters in one product.


- Suitable for single or piecemeal resection
- Versatile
- Reduces need to stock different size snares



Product Code	Opening Width (X1 x Y1 / X2 x Y2)	Diameter	Length
99-05-20-81-23-2	15 x 30mm / 30 x 55mm	2.3mm	230cm

Rotator Snare

- One handed rotation.
- The Rotator Snare is designed for easy, rotating manoeuvrability.
- Braided wire.



Product Code	Opening Width (X x Y)	Diameter	Length
99-05-20-21-22-5	17 x 32mm	2.3mm	230cm
99-05-20-21-23-5	25 x 50mm	2.3mm	230cm



